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A wellness device for  
enhancing our own protective  
forces and correcting the natural  
electromagnetic radiation  
in natural bio-systems

DeVita AP  
Model Base

Methodological recommendations

## Introduction

You are holding a unique, unparalleled anywhere in the world, electromagnetic «DeVita AP Model Base» device – the device for maintenance of high factor of health. Its action is based on the emission of electromagnetic waves, which resonate with (self) oscillations of parasites of different types – from viruses to helminthes. Thus, «DeVita AP Model Base» selectively and deliberately influences causal organisms and destroys them, which allows the human body to function in the natural conditions not weighted down by external factors.

The outstanding effect of «DeVita AP Model Base» use is achieved due to the high power of its electromagnetic field to penetrate into the human body. The oscillation frequency is set with high precision that allows to influence the very microorganisms deliberately, without harming the surrounding organs and tissues.

Specially developed care programs allow to maintain a high health factor in conditions of the most diverse external and internal factors' influence. Each developed course is aimed at influencing a certain microorganism.

The most important advantage of «DeVita AP Model Base» is availability of automatic, specially developed programs for maintenance of high health factor in conditions of the most diverse external and internal factors' influence. Moreover, the owner of the device has a possibility of device programming, using any available preventive programs with «DelInfo» utilization. That enables everyone to select an individual course of preventive care procedures.

## List of programs

- |                           |                        |
|---------------------------|------------------------|
| 1. Anti-Cold              | 12. No. candida        |
| 2. Anti-Flu               | 13. No. chlamydia      |
| 3. Anti-Virus             | 14. No. fungi          |
| 4. Digestive tract health | 15. No. flu, general   |
| 5. Drainage Anti-Tox      | 16. No. helicobacter   |
| 6. Drainage, general      | 17. No. hepatitis A    |
| 7. Help against poisoning | 18. Healthy lymph      |
| 8. Intimate health        | 19. No. herpes         |
| 9. Mexican detoxification | 20. Without toxins     |
| 10. No. allergies         | 21. No. herpes simplex |
| 11. No. bacteria          | 22. No. herpes zoster  |

- |                                    |                                                  |
|------------------------------------|--------------------------------------------------|
| 23. No. mould & fungi              | 27. No. protozoans                               |
| 24. No. nematodes & ascarides      | 28. No. schistosomes (blood-flukes) & trematodes |
| 25. No. parasites, Digestive tract | 29. No. staphylococcus                           |
| 26. No. parasites, general         | 30. No. streptococcus                            |

**Caution!** The device has contraindication to its use by women in the first trimester of pregnancy and by people, who undergone an organ transplantation. In addition, it is not recommended to use the device within the first two months after suffering a myocardial infarction. In case of suffering from serious heart disease, the device must be positioned no closer than 0,5 m from the body.

## 1. Anti-Cold

The causative agents of colds are the following: respiratory viruses, affecting the upper respiratory tract (nose, throat, trachea, larynx). Diseases caused by this group of viruses are airborne and cannot be treated with antibiotics. Complications from a cold can be very dangerous: rhinitis, sinusitis, otitis, etc. Often, when someone catches a cold, bacterial infection can accrue to viral infection due to his (her) low immunity. Therefore, it is advisable to use this program in conjunction with «No staphylococcus» and «No streptococcus» programs. The aim of the program is to reduce inflammation of the mucous membranes of the nose, throat and bronchi, and to activate the local immunity.

**Use:** At the first sign of a cold, during the first 3 days – up to 2-3 times a day. Further – once (1 time) a day for 7-14 days.

## 2. Anti-Flu

Flu is a severe viral infection that affects people of any age. Every year, up to 500 million people become sick with flu. A severe bout of flu may lead to irreversible lesions of the cardiovascular system, respiratory system, central nervous system. The virus is airborne and passed by way of sneezing, coughing and talking. The most common complication of flu is pneumonia. The aim of the program is to destroy different strains of influenza, to remove the skin manifestations.

**Use:** For the prevention of influenza at the epidemic period, the program is used once (1 time) in 2 days. If necessary, once a day, up to 2-3 times a day, alternating with the «No flu, general» program. After that, you should use the program «Drainage, general». In addition, it is possible to use the following programs: «No staphylococcus», «No streptococcus», «No candida», in order to prevent bacterial and fungal contagion; the «Anti-Virus», «Anti-Cold» programs can also be added to this complex.

## 3. Anti- Virus

Viruses represent a transitional form between living and inanimate nature. Since the viruses have genetic material but no cell, they are intracellular parasites. This adds to

the complexity of their destruction: they are under the protection of the cells in which they live. The virus is airborne, household and sexually transmitted. In particular, people with a weakened immune system can be easily infected.

The aim of the program is to fight certain viruses, including cytomegalovirus and Epstein-Barr's virus. Also, the program is used to strengthen the body's immunity and to destroy the virus.

**Use:** On a daily basis over 14 days, during the period of various epidemics, use every other day to enhance the immunity.

#### 4. Digestive tract health

Intestinal infection occurs by ingestion of certain pathogenic viruses or bacteria. A similar infection is usually faced by people of any age, but more often by children and those who take antibiotics. This condition is accompanied by nausea, vomiting, diarrhoea, abdominal pain, and fever. In viral infections other symptoms may also be experienced: a runny nose and redness of the throat, pain at swallowing, lack of appetite, lack of energy.

The aim of the program is to destroy intestinal flu virus (rotavirus) and pathogenic strains of E.coli. The program will help to normalize the situation and to stop dehydration, which is destructive to the body.

**Use:** On the first day use 4-5 times, one after another, in order to achieve sustainable results. Further, use every day: 1-2 times a day for 3-5 days. Alternate with «Drainage, general» and «Anti-Virus» programs. Be sure to drink plenty of liquids.

#### 5. Drainage Anti-Tox

Intercellular matrix is the main place of accumulation of toxins that lead to metabolic disorders, the development of chronic diseases, premature aging of the organism. Drainage means detoxifying tissue types; this helps to remove heavy metals and radionuclides and improves the metabolic processes and the general state of the organism.

**Use:** At severe intoxication – use every day and also in conjunction with «Drainage, general» program.

#### 6. Drainage, general

It is one of the basic programs to clean the body of toxins. Removal of toxic substances from the tissues and organs is performed by the lymphatic system. In the presence of various diseases, a number of harmful substances may increase many times. The lymphatic system has no time to remove them and intoxication can occur again. In addition, illness or injury may affect the process of natural outflow of the wound cavity or organ contents. Excessive accumulation of fluid containing pathogenic organisms may cause severe inflammation and deterioration of the state of a human body.

The program protects the body from unnecessary toxic load, neutralizes and removes toxins from the body, while not impacting the lymphatic system.

**Use:** Once (1 time), after the use of any of the basic anti-parasitic «DeVita AP Model Base» programs.

## 7. Help against poisoning

When enterobacteria or enterovirus enter the organism, they cause an infection to develop, which is popularly called «the poisoning». This infection may be of a hidden nature or possibly more obvious in its manifestation, for example, a headache, muscle spasm in the abdominal area, fever, vomiting, bowel disorders. The consequences may lead to serious injury of the body's organs.

The program helps the body to fight the infection caused by enteroviruses and certain types of Salmonella; it accelerates the recovery process and reduces the probability of complications.

**Use:** The first day 4-5 times (one after another) in order to achieve sustainable results. Further, use it 1-2 times a day, for 3-5 days. Alternate with the programs: «Drainage, general» and «Anti-Virus». Be sure to drink plenty of liquids.

## 8. Intimate health

Genitourinary infections may occur in acute or chronic form.

Chronic infection can eventually cause a narrowing of the urethra, inflammation of the joints, cause infertility and increase the risk of ectopic pregnancy. Periodically, deterioration can occur during a cold or an emotional overload. Infections are mainly transmitted through sexual contact and some of them may be transmitted at birth. Their risk is that they may be asymptomatic in nature.

The program affects the micro-organisms, such as Ureaplasma, Mycoplasma, Toxoplasma.

**Use:** Every day for a month

## 9. Mexican detoxification

Mexican detoxification is the system of complete deep body cleansing, proposed by R. Rife. The program helps to remove metabolic products and toxins from an organism, to improve the drainage function of the lymphatic system, to increase detoxification capabilities of the liver, to activate excretory functions of kidneys and to improve the immune protection.

During the drainage process, the toxins transfer from the matrix through the lymphatic system into the bloodstream. From the bloodstream the toxins are removed by the organs of the excretory system, mainly by the liver and kidneys. In the course of drainage, there is a risk of overload to the organs of the excretory system, which may develop a dysfunction of these vital organs. For this reason the «Mexican detoxification» program includes frequencies, which support the organs of the excretory system.

**Use:** Every day, during one month, at the absence of serious diseases of the liver, kidneys, and lymphatic system. In case of having any of those, it is recommended to use the «Complete cleansing» automatic program for genital cleansing of the body.

## 10. No allergies

An allergy is a hypersensitivity of the organism to any substance, food, wool, dust, pollen, chemical and infectious agent. Almost half of the world's population suffers from some sort of allergy. Often the reason for the allergy is an improper diet and unhealthy lifestyle.

The program relieves the basic symptoms of the allergy, alleviates the human condition, and helps to stabilize the immune system.

**Use:** Every day – when it is indicated, and to ease the symptoms. Use 2-3 times a day, if particularly exacerbated symptoms persist.

## 11. No bacteria

The most dangerous in the development of diseases caused by a bacterial infection is a bacterial toxin excretion (poisons), which creates an inflammatory action, intoxication and damages internal organs.

The program is used in the presence of infecting by bacteria of *Klebsiella*, *Listeria*, and *Pneumococcus* to reduce the severity of symptoms of intoxication and to facilitate the human condition, at the point of a bacterial infection of the digestive tract or of the respiratory system.

**Use:** Every day, until signs of the infection have disappeared.

## 12. No candida

*Candida* is part of the normal microflora of intestines, vagina, or mouth. But with lowered immunity *Candida* begins to proliferate and harm the body. There is also pathogenic *Candida*, which enters the organism with food, through household items, after the use of potent antibiotics, after sexual contact with an infected partner. Symptoms of infection are manifested in the form of allergic rashes on the skin, thrush, pain and burning in the genital organs, dysbacteriosis.

The program allows the preventing of candidiasis.

**Use:** 1-3 times a day for 2-4 weeks. Use in the complex with «No protozoans» program.

## 13. No chlamydia

Chlamydia causes infectious sexually transmitted diseases. Chlamydia is often asymptomatic or with minimal symptoms. It affects mainly urogenital organs. Women are more susceptible to this infection; children can become infected when passing through the mother's birth ways. Chlamydia provokes serious complications, including infertility, lesions of the urogenital organs, joints and eyes if the infection is not treated in time.

**Use:** Every day, for a period of 2-3 months.

## 14. No fungi

The main danger of fungal disease is that areas affected by fungus open the way for the pathogens of various infections. The most common fungal infections affect the skin and its derivatives (hair, nails), and mucous membrane of the internal organs. Treatment requires an integrated approach and good hygiene. The program actively fights the fungi: such as Trichophyton, Mikrosporum, Mucor.

**Use:** In order to prevent invasion by fungal infections, it is recommended to use the program each time after visiting the sauna, pool, beach, after contact with an infected person. If there is a fungal infection in the body: use 1-2 times a day during a 1 month period.

## 15. No flu, general

The flu is a disease that affects the upper and lower respiratory tract; it is accompanied by severe intoxication and can lead to serious complications. Flu begins acutely and quickly (48 hours from the moment of infection until the first signs). The main symptom of intoxication is a raised body temperature (40° C and above). This temperature can last for more than 5 days which causes shortness of breath or respiratory failure, impairment of consciousness, seizures. For people with weakened health, the elderly, small children, this temperature leads to an infectious-toxic shock (a condition when the vital organs are denied. By using «No flu, general» program, you will be able to alleviate the condition of the patient and by complex use of several programs against influenza viruses you can quickly restore health and prevent the complications.

**Use:** Use at the first sign of flu symptoms with increased body temperature. Use the program 1-3 times (one after another), or alternate it with «Anti-Flu» program. Be sure to drink plenty of liquids.

## 16. No helicobacter

Helicobacter Pylori is a microorganism, living in the stomach and duodenum of humans and provokes the development of dangerous diseases: gastritis, ulcers, hepatitis, etc. According to WHO (World Health Organisation), 60% of the population is infected with Helicobacter Pylori. Infection occurs easily: when a family member is infected, the chance of other members to be infected is 95%. Transmission of the parasite occurs as a result of the general use of personal care products, tableware and it is also airborne. This microorganism may cause such symptoms as heartburn, bloating, pain, heaviness in the stomach.

**Use:** Every day, for 2 weeks.

## 17. No hepatitis A

Hepatitis A virus is transmitted through dirty hands, food, water, including sea water. It is especially common in Asia and Africa. When planning travel to these countries

you should take appropriate precautions: observe hygiene, drink only bottled water and use «No hepatitis A» program for the prevention and at the first suspicion of the presence of infection in you or in other people.

**Use:** Every day up to 2-3 times a day, for at least a month.

## 18. Healthy lymph

The lymph is the intercellular fluid found in the lymphatic system. The role of the lymph is particularly important as it contributes to the maintenance of water balance, homeostasis of the microenvironment of cells and has the potential to eliminate pathological microorganisms and toxins from the tissues. The inadequate circulation of the lymph nodes results in the non-response of the lymph function, leading to a pathological process.

The program aims to normalize lymphatic drainage, enhance lymphatic vessels, improve lymph nodes, which leads to decreased edema and better release of metabolites from the tissues. The use of this program, apart from the above, increases the defense of the organization.

**Use:** every day for one month.

## 19. No herpes

More than 90% of all people in the world are the carriers of different types of the herpes virus. You can become infected in many ways: through saliva, sex, unhygienic items, as well as airborne. At the beginning the disease may be asymptomatic, then the following characteristic symptoms may appear: headache, malaise, fever, herpetic eruptions. The program will help to reduce the recovery period and reduce inflammation and pain at various herpes infections.

**Use:** 1-2 times a day for a month.

## 20. Without toxins

The program contributes to the removal of toxins, improves drainage (detoxification) functions, resulting in improved immunity and increased vital energy. The program can be used after anti-parasitic programs, and independently for one month. Unlike other detox programs, it does not consume the body's energy resources.

**Use:** every day for one month.

## 21. No herpes simplex

Every third person is a carrier of the herpes simplex virus. This virus is localized in the mucous membranes of the lips, mouth, and genitals. It is transmitted mainly by contact with the fluids of the infected person's body (blood, saliva, semen, mucus secretion). If a person has been infected, herpes continues to exist in his body in a latent form, which exerts in the form of a rash when the immunity is decreased:

during hypothermia or overheating, during infectious diseases or during menstruation in women. The program effectively reduces the herpetic complications and prevents occurrence of the infection signs in the mucous.

**Use:** At the first signs of rash every 3 hours, and then 1-2 times a day for a month.

## 22. No herpes zoster

The disease, caused by herpes zoster virus, is called «shingles». If a person is ill with chicken pox (a form of zoster virus), the virus remains in the body in a latent form, and hides in the nerve cells. When person becomes ill with cold or chronic diseases, the virus enters the skin surface and forms on it the characteristic rash and severe pain on the path of any sensory nerve.

The aim of the program is to destroy the virus and to remove the main signs of its manifestation.

**Use:** Everyday up to 3-4 times a day, until the disappearance of symptoms of the disease, at least during 1 month. Use in complex with «No herpes» program.

**Caution!** The program is not used when a person is ill with varicella.

## 23. No mould & fungi

Modern heated apartments or private houses are the perfect abode for mould and fungi. Favorite habitats of mould are bathrooms, libraries, poorly ventilated premises. Fungi can enter the body through food, water, air. Mould penetrates the skin barrier wherever there are suitable conditions for it: humidity, temperature, organic residues in the form of mud or sebum. Fungi and mould cause diseases of the respiratory organs, digestive tract, skin and allergy.

The program is used to prevent fungal infections and to reduce manifestations of allergic reactions.

**Use:** 1-3 times a day, during a 4 week period. Use in complex with «No fungi» program.

## 24. No nematodes & ascarides

The eggs of helminthes, such as nematodes and ascaris, are ingested with unwashed fruits and vegetables, undercooked meat, through dirty hands. Living in their host organism, helminthes eat their tissues, process their nutrients (that is why there is an unhealthy appetite), and produce toxins. The infected person starts rapidly to lose weight, his/her skin peels, and he/she has hair loss. The most common helminthes are roundworms (ascaris). They can cause damage to the wall of the small intestine, intestinal obstruction, parotitis, bronchitis, focal pneumonia, hemoptysis, diseases of the liver and pancreas. Parasitism of ascarids and hookworms causes an allergic reaction, which manifests in the form of asthma and skin rash. Under the action of the program the helminthes stop eating, which leads to their death and to the natural withdrawal of the parasites from the digestive tract.

**Use:** On a daily basis. In total you should have at least 10 sessions. Repeat the course after 14 days. People who often communicate with animals and small children are advised to carry out preventive measures regularly: once (1 time) every quarter (every 4 months).

## 25. No parasites, Digestive tract

The presence of parasites in the digestive tract and in the liver causes different diseases: such as cholangitis, biliary dyskinesia, cholecystitis, hepatitis. The presence of parasites in the gut causes inflammation of the intestinal wall and reduced absorption of nutrients. Parasites release toxins which cause pain in the joints and muscles, affect the central nervous system, contribute to the emergence of depressive state, nervousness, restlessness, and sleep disturbances.

The aim of the program is to cleanse the body, liver and pancreas from the gastrointestinal parasites.

**Use:** The program is used on a daily basis, during 7-10 days. Repeat the course after 14 days.

## 26. No parasites, general

Wellness of the body without decontamination of parasites does not give a positive effect. The parasites are a trigger for many chronic diseases: cholecystitis, cholelithiasis, pancreatitis, colitis, diabetes, asthma, dermatitis. Symptoms such as chronic fatigue, irritability and anxiety, hyperactivity in children, anemia, brittle nails and hair, problem skin, headaches, eating disorders, low immunity, can be signals of this parasitosis. All parasitic diseases are dangerous because they are perfectly disguised as other diseases.

**Use:** The program is used on a daily basis. In total you should have at least 10 sessions. Repeat the course after 14 days.

## 27. No protozoans

Protozoa live in the body of every fifth person; infection occurs through dirty hands or food and water. The infected person feels severe abdominal pain, bloating, constant feeling of fatigue, depression. His liver and pancreas can suffer.

The program affects the most common types of protozoa: such as giardiasis, amoebas.

**Use:** 1-2 times a day, for 10-14 days. If re-infection occurs, it is necessary to identify the source of infection.

## 28. No schistosomes (blood-flukes) & trematodes

Trematodes are flat worms which cause serious infectious diseases. Trematodes can parasitize in the blood (schistosoma), liver, intestines, and lungs. The presence

of these parasites in the body can also become the cause of severe allergies in the form of itching, biliary dyskinesia, increased liver size, headaches and abdominal pain. Infection occurs by eating poorly washed greens or lettuce, not fully cooked meat of cattle, fish.

The aim of the program is to remove most trematode species from the organism (including schistosoma), irrespective of their stage of development and their location in the body.

**Use:** On a daily basis, for 14 days.

## 29. No staphylococcus

According to statistics, 20% of the population are constant carriers of staphylococcus. The bacteria live on the skin and in the upper part of the airways. Staphylococcus aureus is a frequent pathogen of infectious diseases, which affects the skin, soft tissue, upper respiratory tract, musculoskeletal system. Staphylococcus causes complications in the presence of flu, festering wounds, food poisoning.

The program works for common types of staphylococci, and is effective in the treatment of diseases caused by these bacteria.

**Use:** On a daily basis, up to 2 weeks, depending on the manifestations of infection.

## 30. No streptococcus

Streptococcus is a widespread group of microorganisms, including pathogens. They live on human skin, in the respiratory tract, in the intestine. They cause acute inflammation of the upper respiratory tract, urinary tract, abdominal organs, skin and soft tissues. Streptococcal infection is more severe than staphylococcus and can lead to disease of the joints, heart, kidneys, upper respiratory tract and others. The program is effective in the presence of angina, otitis, sinusitis and can be used both in complex therapy and in monotherapy.

**Use:** On a daily basis, up to 2 weeks, depending on the manifestations of infection.

## Automatic programs

| No in device | Name of the program   | Aim of program                                                                                                                                                                                                   | Use                                       |
|--------------|-----------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------|
| 31           | Mode 1 - automatic    | Oppresses the activity of viruses, inhibits the vital activity of protozoal microorganisms. Optimal prevention of diseases when changing climatic zones, like visiting Asian countries.                          | Daily, during 1st week of the first month |
| 32           | Mode 2 - automatic    | Removes toxic load caused by strepto- and staphylococcus. It is effective both for colds and for the prevention of coccal infections.                                                                            | Daily, during 2nd week of the first month |
| 33           | Mode 3 - automatic    | Promotes the removal of bacterial complications of the gastrointestinal tract and urogenital system.                                                                                                             | Daily, during 3rd week of the first month |
| 34           | Mode 4 - automatic    | It displays mycotoxins formed by molds and yeast fungi. Reduces the risk of accumulation of carcinogenic factors.                                                                                                | Daily, during 4th week of the first month |
| 35           | Cleansing - automatic | It stimulates the removal of toxins and the normalization of lymphatic drainage, strengthening the lymphatic vessels. Helps removing the excess fluid from the intercellular space, it has an antiseptic effect. | Daily, during the second month            |

**Application:** in order to avoid a large toxic load, it is not recommended to use the «Cleaning» program and the «Complete Cleansing» automatic program of DeVita Ritm in one month.

**Attention!** If during the use of the programs of the first month nausea, headache, perspiration, or other symptoms of intoxication are observed, then additional application of this or other mode has to be carried out within the next 3 months. Each of them is involved in the course taking into account the life cycle of the infectious organism.

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Deta-Elis Europa GmbH