DEVITA COSMO

DIGITAL DEVICE FOR HEALTH IMPROVING

BEAUTY AND YOUTH

This is a unique device to conduct cosmetic procedures by your own. Revolutionary approach to solving skin problems and a person exterior. The device ensures high efficiency due to simultaneous action of light quantum on nerve terminals and electromagnetic revitalizing action on a body organs and systems.

Beauty and rejuvenation impact is based upon bioresonance method and reverse-reflex action of quantum radiation.

EFFECTIVE. HARMLESS. PHYSIOLOGICAL. EASY-TO-USE.

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DIGITAL DEVICE FOR HEALTH IMPROVING

BEAUTY AND YOUTH

Takes care of skin, hair and nails beauty

Controls performance of organs and systems

Improves skin exterior

Rehabilitates endocrine system of men and women

Triggers rejuvenation

Eliminates inflammations and age spots

Favors rhytides effacement

Reactivates skin tonus and elasticity

DEVITA COSMO

DIGITAL DEVICE FOR HEALTH IMPROVING

BEAUTY AND YOUTH

20 PROGRAMS

QUANTUM AND
ELECTROMAGNETIC
TECHNOLOGY SYNERGY

POSSIBILITY TO USE
WITHOUT ACCESSORY


DEVITA COSMO

DIGITAL DEVICE FOR HEALTH IMPROVING

LIST OF PROGRAMMS


1. Wrinkle-free skin
2. Skin tonicity
3. Female health
4. Male health
5. Liver performance
6. Live peacefully
7. Ventricle regulation
8. Digestion
9. Like clockwork
10. Nephros activation
11. Thyroid body
12. Heart and vessels
13. Spleen performance
14. Key to life
15. Bright eyes
16. Clear skin
17. Beautiful braid
18. Nail artist
19. Flaws-free face
20. Soft heels

QUANTUM ACCESSORY DEVITA COSMO

DIGITAL DEVICE FOR HEALTH IMPROVING

PRIMARIES OF LIGHT


DEVITA COSMO OPERATING PRINCIPLE

DIGITAL DEVICE FOR HEALTH IMPROVING


EVERY HEALTHY ORGAN
HAS ITS OWN
REFERENCE
ELECTROMAGNETIC
OSCILLATION
FREQUENCY

OSCILLATIONS OF A
FRAGILE ORGAN
DEVIATE FROM
REFERENCE
OSCILLATION
FREQUENCY OF A
HEALTHY ORGAN

DEVITA COSMO OSCILLATES REFERENCE ELECTROMAGNETIC FREQUENCIES OF HEALTHY ORGANS AND TUNES A BODY TO A HEALTHY FREQUENCY WITH THE USE OF RESONANCE EFFECT

Completely harmless Influence on the cells of living organisms including human body cells

DEVITA COSMO

DIGITAL DEVICE FOR HEALTH IMPROVING

QUANTUM ACCESSORY PRINCIPLE OF OPERATION

Light quanta possess wide spectrum distribution which consists of emissions of different wave and frequency ranges.

Every color makes definite impact on the organs or systems through cutaneous receptors in organ projection area.

Such impact is based upon Zakharyin-Ged inverse reflex feedback on facial zones.

Due to inverse impact on nerve terminals and synthesis of biologically active substances enhancement, quantum radiation renews balance of body systems and favors their functions recovery.


BIORESONANCE

DIGITAL DEVICE FOR HEALTH IMPROVING

THE WORD RESONANCE DERIVES FROM THE LATIN WORD *RESONO*
THAT MEANS TO RESPOND

RESONANCE is system oscillation amplitude increase as a result of external impact accompanied by external impact frequency matching with internal oscillation frequency in the system.

BIORESONANCE is an enhancement of biological system feedback (response) to external periodic effect in a narrow frequency range.

PERMITTING DOCUMENTS

DIGITAL DEVICE FOR HEALTH IMPROVING

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CONCLUSION: IN REFER TO SGS TEST REPORT NO. HIF40005
NO DEVIATION FROM THE STANDARDS HAVE
BEEN IDENTIFIED. THEREFORE THE TESTED
PRODUCT IS IN COMPLIANCE WITH THE
APPROPRIATED STANDARDS.

This Verification of Compliance has been prepared on the basis of the results of analysis, performed by Laboratory of SGS GERMANY OMBH, Munich, Germany. It is an integral part of the compliance certificate and is valid only in the presence of the relevant quality standard. The CE marking on devices shall only be affixed, under the responsibility of the manufacturer, after compliance with CE Declaration of Conformity and verification with relevant CE Standards.

Month: June 21, 2016
City: Detmold

SGS Germany GmbH
Munich, Germany

SGS

CONFIRMATION OF COMPLIANCE

CUSTOMER: DETA-ELIS EUROPA OMBH
JUSTUS-LIEBIG-STR. 2-4
3609 KUNZEL
GERMANY

SERVICE SUPPLIER: SGS GERMANY OMBH
EMC LABORATORY MUNICH
HOFMANNSTRASSE 59
81739 MUNICH
GERMANY

MODELS: DeVita AP Model Mini; DeVita RITM Model Mini;
DeVita Energy

SPECIFICATIONS: EN 55014-1:2006 + A1:2009 + A2:2011
EN 55014-2:1997 + A1:2001 + A2:2008
EN 61000-3-2:2006 + A1:2009 + A2:2009
EN 61000-3-3:2008

CONCLUSION: IN REFER TO SGS TEST REPORT NO. HIF40005
NO DEVIATION FROM THE STANDARDS HAVE
BEEN IDENTIFIED. THEREFORE THE TESTED
PRODUCT IS IN COMPLIANCE WITH THE
APPROPRIATED STANDARDS.

This Verification of Compliance has been prepared on the basis of the results of analysis, performed by Laboratory of SGS GERMANY OMBH, Munich, Germany. It is an integral part of the compliance certificate and is valid only in the presence of the relevant quality standard. The CE marking on devices shall only be affixed, under the responsibility of the manufacturer, after compliance with CE Declaration of Conformity and verification with relevant CE Standards.

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